


Button Injection Site



- Large & Arc-surface design prevents the bacterial buildup and reduces alcohol retention when swabbing.
- Transparent device provides immediate detection of blood flashback.
- Non-removable compressed port instantly reseals, permitting repetitive injections with relatively large gauge needles.
- Synthetic Polyisoprene Site to eliminate latex sensitivity.



1)
(

Part Number	Characteristic / Material	Description
SA0137	Toray ABS T-900 Site- Latex Free Polyisoprene	Button Injection Site I.D. 2.08mm/O.D. 3.2mm
SA0138	BASF TERLUX ABS 2802 Site- Latex Free Polyisoprene	Button Injection Site I.D. 2.08mm/O.D. 3.2mm